Instructions for Use

Dear customer, thank you for purchasing an EMO product. By choosing an orthosis from the NeoActiv line you have purchased a quality product which meets high medical-health standards. Please read these instructions carefully. If you have any questions please contact your doctor or your specialised establishment.

Regulations

The medical device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices and the respective national regulations. Class I medical device, non-sterile, without measurement function.

Description and Characteristics

The NeoActiv orthoses are made of medical/health quality neoprene, a thermal insulator with elasticity and viscosity indexes selected for therapeutic use. This material conserves and increases the body's temperature in the area it is fitted, speeding up and shortening the heating time, increasing the elasticity of the muscles and tendons, increasing the effectiveness of the rehabilitation, providing therapeutic heat and preventing injuries. Its even elasticity improves its adaptation and makes it easier to fit it. Its excellent compression does not limit mobility. It relieves pain and accelerates the recovery process at the same time as it reduces the oedema.

Indications

The indications listed below are the most important. All the orthoses are indicated for thermocompression, prevention, protection, stabilisation, reduction of oedemas, inflamatory processes, post-traumatic and post-surgical rehabilitation and protection during sporting activities. Wrist supports: slight articular trauma-tisms, tendonitis, rheumatic pain, rhizar-throsis, bursitis, sprain. If a splint is incorporated, immobilisation and stabilisation of the joint. Elbow supports: muscular and ligament injuries (bursitis, fasciitis, rheumatic pain), epicondylitis and medial epicondylitis. Shoulder supports: subluxation of the scapular-humeral joint, instabilities, painful processes. Lumbar supports: lumbagos, muscular hipotony, postural and post-surgical rehabilitation and activities which

require a continuous mechanical overload of the lumbar-sacral zone. Trousers: overloading and inguinal adductor injuries. High supports: torn muscles. Leg supports: torn muscles and protection of the tibial crest. Knee supports: arthritis of the knee, tendonitis, inflamatory processes, slight ligament strains, patellar centring, injury of the patellar tendon. Knee supports with side bars: instability, arthritic and arthrosic pains and contusions. Knee supports with articulations: chronic instability, control of flexion-extension, post-traumatic and post-surgical rehabilitation. Ankle supports: tendonitis, postural control, mid-lateral instability (sprains), avoiding relapses whilst playing sport.

Instructions for Use and Fitting

The product must be adapted by the Orthopaedic Technician or by the doctor. Wrist supports: they are fastened and adjusted using Velcro. If a splint is incorporated it must be adapted to the patient's anatomy. Elbow supports: if straps are incorporated for epicondylitis fit when the muscles are relaxed and then adjust. Knee supports: the patellar opening must be correctly centred on the patella. If bars or articulations are incorporated, these must be parallel on both sides of the joint. Then adjust the Velcro straps. Infrapatellar strap: fit between the patella and the lateral condyle of tibia. Adjust symmetrically using the pull strips. Ankle supports: fit the ankle support with the ankle in neutral position of 90°. Then, if the pronator straps are incorporated, adjust them for the dressing in an 8-shape. If it incorporates mid-lateral plates, place them in the side pockets.

Precautions

Follow the manufacturer's directions for proper orthosis placement. Follow the instructions or recommendations of the prescriber or healthcare personnel. For a personalized adaptation, consult a specialist. If there is ongoing pain, it is advisable to contact a healthcare professional. Fit ensuring that the compression is not excessive. If in any doubt when choosing the size always opt for the larger size. Only use during periods of activity. Do not expose the item to naked

flames or heat sources at very high temperatures. Guarantee within the framework of the legal regulations.

Washing and Maintenance

Wash by hand with warm water and neutral detergent. Whilst washing, if the item contains Velcro strips, they must be stuck together or kept closed. Dry at room temperature. Do not hang out or iron the item. Contact with ointments, creams, alcohols and solvents can deteriorate the item.

Quality Certificates

All the materials used have been tested and comply with applicable European regulations, and are high quality and offer the maximum safety and comfort during use. Öko-Text Standard 100 Certified: Organic textile label awarded after the drawing up a catalogue of analysis criteria which guarantees the absence of dangerous substances for the health. Certified as the absence of harmful substances in all the materials with which the item is made. Tested by the Institut Fresenius: Chemische und Biologische Laboratorien: According to the International Contact Dermatitis Research Group criteria the neoprene used in NeoActiv achieves excellent tolerance results in all kinds of patients, even in those people with very sensitive skin, to whom it does not cause any type of skin reaction.

Notice

To dispose of the product, use the solid waste container. For the disposal of the container, depending on its material, use the paper / cardboard container or the plastic container. In any case, you must strictly comply with the regulations of the country. Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and / or the patient are established.

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