NightSplint®

EN Splinting System
DE Lagerungsschiene
ES Férula estabilizadora
FR Attelle de positionnement
IT Ortési da riposo
NL Immobilisatiebrace
DA Lejringsskinne
SV Viloskena
PL Szyna unieruchamiająca
SL Položajna opornica
HR Imobilizacijska ortoza
PT Tala de posicionamento

DARCO

NightSplint® | Splinting System

Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This aid assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product with injured, irritated, damaged or broken skin. You should use bandages or socks, for example.

This aid is intended for single-patient use only.

Purpose

The NightSplint® splinting system is designed to stabilize dorsal extension which consequently results in a slight stretching of the Achilles tendon, the calf muscles as well as the plantar fascia.

Indications / Contraindications

Indications

- > Shortening of the plantar fascia
- > Plantar fasciitis
- > Achilles tendonitis
- > Prophylaxis of talipes equines
- > Postoperative immobilization

Contraindications

- > Plantar flexion contracture exceeds 10°
- > No control with severe plantar flexion spasticity

Design features

- > Two heel wedges for adjusting dorsal flexion
- > Three hook and loop fastener straps for selective customization
- > Smooth upper material – high comfort
- > Non-slip sole

Technical data

Sizes	Shoe size	Length	Height	Width
S	35 – 36	23.0 cm	33.0 cm	12.0 cm
M	37 – 39	23.5 cm	35.0 cm	13.0 cm
L	40 – 43	28.0 cm	37.0 cm	14.0 cm
XL	44 – 49	34.0 cm	38.0 cm	16.5 cm

Universal fit for right and left fit.

Colour	Black
Material	Splinting System: Nylon Liner: Polyester foam 2 Heel wedges: EVA

Positioning the foot into the splinting system:

1. Open the hook and loop fastener straps at the front in the toe area underneath the splint. Remove the two heel wedges from the wedge pocket.
2. Enter the foot into the splinting system with the knee bended in order to realize the angle of maximum dorsiflexion.
3. Fix the strap across the instep with the hook and loop fastener strap in order to optimize the position of the heel.
4. Stretch the knee in order to check the fitting of the strap across the instep (If the heel is slipping inside the splint the strap is too loose, if the calf is squeezed in, the strap is too tight. Readjust and check again).
5. After checking the correct fit of the strap across the instep, the other two hook and loop fastener straps can be closed across ankle and shinbone in order to ensure a secure and comfortable fit.

Inserting the heel wedges into the splinting system:

1. Insert the heel wedges only under supervision of a physician or healthcare specialist.
2. With the foot in the splint open the strap located underneath the splint and initially insert one wedge thin end first.
3. Try and test this stabilizing orthosis for one hour. In case of skin irritation, changes in color or increased pain discontinue use and contact your physician or healthcare specialist immediately.
4. If during the night no pain or discomfort occurs and your medical healthcare specialist agrees, the second wedge may be inserted the same way as the other.

Combination with other products

No specifications

Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling / cleaning. The medical device is designed to be used for the duration of use prescribed by a doctor without showing a reduction in safety or performance. After the period of use prescribed by your doctor, check the function of the product together with your doctor or medical provider.

Maintenance

The product does not require maintenance.

Cleaning instructions

1. Remove the heel wedges.
2. Open the hook and loop fastener straps in the heel area of the product.
3. Take out the plastic body by bending together the two shell parts.
4. For cleaning use mild detergents only. The liner is only to be hand washed and air dried.

Disposal

Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. Dispose of this product after use in accordance with applicable local laws and regulations.

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DARCO International, Inc.
ALBO-Healthcare GmbH