# DAW

## (DARCO Achilles Wedges)

- EN Heel Wedges for Walkers and Foot-Positioning Orthosis
- DE Fersenkeile für Walker und Fußlagerungsorthese
- ES Cuñas de talón para Walker y Ortesis de posicionamiento del pie
- FR Cales de talon pour Walker et Orthèse de positionnement du pied
- IT Zeppa per tallone per Walker e Ortesi di posizionamento del piede
- NL Hielwiggen voor Walkers en Voetimmobilisatie-orthese
- DA Hælkiler til Walker og Fodlejningsortose
- SV Hälkilar för Walker och Fotavlastningsortos
- PL Podpiętki do ortozy typu Walker oraz Orteza pozycjonująca stopę
- SL Klinasti poviški za peto za pripomoček Walker  
in Ortoza za imobilizacijo stopala
- HR Ulošci za petu za ortoza za Walker  
i Ortoza za imobilizaciju stopala
- PT Cunhas de calcanhar para Walker  
e Órtese de posicionamento do pé

# DARCO

## EN Instructions for use

### DAW (DARCO Achilles Wedges)

Heel Wedges for Walkers and Foot-Positioning Orthosis

#### Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your doctor, specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, they will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation, etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This device assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility.

Do not wear the product in direct contact with injured, irritated, damaged or broken skin. You should use bandages, dressings or socks. You must not drive a vehicle while you are wearing this product.

Use extreme caution when climbing stairs (e.g. also stepping up onto curbs) and walking on uneven surfaces. Be particularly cautious when walking in wet, slippery and snowy conditions. Walk more slowly than usual.

This device is allowed for single-patient use only. It does not need to be cleaned or reprocessed between uses.

## Purpose

Heel wedges for defined adjustable positioning of a lower leg / foot orthosis or foot-Positioning Orthosis.

## Indications / Contraindications

### Indications

- > Achilles tendon ruptures
- > Postoperative use for Achilles tendon suture
- > Postoperative use for removal of Haglund's exostosis
- > Tendinitis in the Achilles tendon area
- > Plantar flexion contracture of more than 10°
- > Shortening of the calf muscles

### Contraindications

- > Forefoot conditions
- > Plantar foot ulcers

## Design features

DAW heel wedges consist of 3 wedges in different lengths. The wedges are placed on top of each other according to the length of the foot and the desired angle of plantar flexion and fixed to the insole in the heel area of the walker or positioning orthosis. The wedges come in one standard size. They can therefore be used universally with all sizes of walkers or positioning orthosis. They can be used on the right and left.

## EN Use

1. Open the walker /foot-positioning orthosis and leave the insole inside the walker /foot-positioning orthosis .
2. Place the wedges in the heel area of the walker /foot-positioning orthosis, corresponding to the degree of plantar flexion that is therapeutically required. The wedges are placed on the insole of the walker /foot-positioning orthosis.
3. A wedge has a height of approx. 11 mm. Depending on the number of wedges used, the height varies and the angle of plantar flexion can, thus, be adjusted individually. 3 wedges are equivalent to approx. 120° slant.
4. Fix the wedges to the insole of the walker /foot-positioning orthosis using the adhesive dots.
5. Close the walker /foot-positioning orthosis .
6. Depending on the course of therapy, the wedges can be removed and, thus, the height adjusted.

## Combination with other products

Relief Insert® 2.0  
Relief Insert® Walker 2.0  
Relief Insert® Walker 2.0 Achill

## Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling / cleaning.

The medical device is designed to be used for the usual duration of use prescribed by a doctor without showing a reduction in safety or

performance.

After the period of use prescribed by your doctor, check the function of the product together with your doctor or specialist dealer.

## Technical data

**Size:** Universal  
**Contents:** 3 wedges  
**High wedge:** Approx. 11 mm per wedge; equivalent to approx. 120° plantar flexion  
**Colour:** Black  
**Material:** EVA

## Maintenance

The product does not require maintenance.

## Care instruction

If the wedges are slightly soiled, they can be wiped with a damp cloth. Do not dry the wedges in a dryer or on a radiator. Leave the wedges to dry at room temperature.

## Disposal

 Dispose of after use.  
Please dispose of the packaging in an environmentally friendly manner.  
Dispose of this product after use in accordance with applicable local laws and regulations.

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ALBO-Healthcare GmbH