WCS[®] Light

EN Closed Wound Care Shoe
DE Geschlossener Verbandschuh
ES Zapato postquirúrgico cerrado
FR Chaussure fermée pour le soin des plaies
IT Calzatura ortopedica terapeutica chiusa
NL Gesloten verbandschoen
DA Lukket bandage- og sårsko
SV Sluten sårbehandlingssko
PL Zamknięty but do terapii ran
SL Zaprt čevelj za nego ran
HR Zatvorena ortopedska cipela
PT Sapato fechado para terapia de feridas

DARCO

1

2


3

4


5

6

7




EN Date of manufacture
DE Herstellungsdatum
ES Fecha de fabricación
FR Date de fabrication
IT Data di fabbricazione
NL Productiedatum
DA Fremstillingsdato
SV Tillverkningsdatum
PL Data produkcji
SL Datum izdelave
HR Datum proizvodnje
PT Data de fabrico

EN Manufacturer
DE Hersteller
ES Fabricante
FR Fabricant
IT Fabricante
NL Fabrikant
DA Producent
SV Tillverkare
PL Producent
SL Proizvajalec
HR Proizvođač
PT Fabricante

EN Swiss Authorised Representative
DE Schweizer Bevollmächtigter
ES Representante autorizado en Suiza
FR Représentant autorisé en Suisse
IT Rappresentante Autorizzato Svizzero
NL Gevolmachtigde vertegenwoordiger Zwitserland
DA Autoriseret repræsentant i Schweiz
SV Auktoriserad representant i Schweiz
PL Autoryzowany przedstawiciel na terenie Szwajcarii
SL Pooblaščen predstavnik v Švici
HR Ovlašteni zastupnik u Švicarskoj
PT Representante autorizado suíço

A**B****C****D****E****F****CE**

- EN Certification mark declaring this product's compliance with EU legislation
 DE Zertifizierungszeichen, das die Übereinstimmung dieses Produktes mit der EU-Gesetzgebung bestätigt
 ES Marca de certificación que declara la conformidad de este producto con la legislación de la UE
 FR Marque de certification attestant de la conformité de ce produit à la législation de l'UE
 IT Marchio di certificazione che dichiara la conformità del prodotto con la legislazione UE
 NL Certificeringskenmerk dat aangeeft dat dit product voldoet aan de wetgeving van de EU
 DA Certificeringsmærke, der erklærer dette produkts overensstemmelse med EU-lovgivningen
 SV Märket deklarerar att produkten uppfyller lagarna i EU
 PL Znak certyfikacyjny deklarujący zgodność tego produktu z przepisami UE
 SL Certifikacijska oznaka, ki potrjuje skladnost tega izdelka z zakonodajo EU
 HR Certifikacijska oznaka koja označava sukladnost ovog proizvoda sa zakonodavstvom EU
 PT Marca de certificação que declara a conformidade deste produto com a legislação da UE

**UK
CA**

- EN Certification mark declaring this product's compliance with UK legislation
 DE Zertifizierungszeichen, das die Übereinstimmung dieses Produktes mit der Gesetzgebung des Vereinigten Königreiches bestätigt
 ES Marca de certificación que declara la conformidad de este producto con la legislación de GB
 FR Marque de certification attestant de la conformité de ce produit à la législation de l'RU
 IT Marchio di certificazione che dichiara la conformità del prodotto con la legislazione del Regno Unito
 NL Certificeringskenmerk dat aangeeft dat dit product voldoet aan de wetgeving van de VK
 DA Certificeringsmærke, der erklærer dette produkts overensstemmelse med UK-lovgivningen
 SV Märket deklarerar att produkten uppfyller lagarna i Storbritannien
 PL Znak certyfikacyjny deklarujący zgodność tego produktu z przepisami Wielkiej Brytanii
 SL Certifikacijska oznaka, ki potrjuje skladnost tega izdelka z zakonodajo Združenega kraljestva
 HR Certifikacijska oznaka koja označava sukladnost ovog proizvoda sa zakonodavstvom UK
 PT Marca de certificação que declara a conformidade deste produto com a legislação do Reino Unido

EN Instructions for use

WCS® Light | Closed Wound Care Shoe

Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your doctor, specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation, etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This device assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product in direct

contact with injured, irritated, damaged or broken skin. This device is not recommended for direct skin contact. You should use bandages, dressings or socks.

You must not drive a vehicle while you are wearing this product. Use extreme caution when climbing stairs (e.g. also stepping up onto curbs) and walking on uneven surfaces. Be particularly cautious when walking in wet, slippery and snowy conditions. Walk more slowly than usual.

Note: This device is allowed for single-patient use only. It does not need to be cleaned or reprocessed between uses.

Purpose

The Wound Care Shoe is mainly used for minor ulcerations. This enables specific offloading and pressure redistribution of the foot.

Indications / Contraindications

Indications

- > For treatment of diabetic foot syndrome
- > For healing of open wounds and minor ulcerations
- > For pressure relief
- > During the postoperative phase after toe and/or forefoot surgery

Contraindications

- > No contraindications are known.

Design features

The flexible upper material optimally protects the foot from external influences and prevents dorsal pressure points. The special sole construction achieves optimum pressure redistribution. The multi-layer insole can be adapted individually to the patient to achieve selective offloading. It can be worn as an everyday shoe both indoors and outdoors.

Technical data

Sizes

XS (35 – 36), S (37 – 38),
S / M (39 – 40), M (41 – 42), L (43 – 44),
XL (45 – 46), XXL (47 – 48)

Available in pairs.

Material

Outsole: Rubber, EVA, bamboo
Insoles: Plastazote®, EVA, PU
Upper material: Polyester, Plastazote®,
synthetic leather, PU

Colour: Black

Combination with other products

WCS®/WCS® Light Insoles

Maintenance

The product does not require maintenance.

Modifications may only be carried out by specialised medical personnel. Determining the correct shoe size is very important for the functionality of the shoe.

1. Determine the correct shoe size from the insole length.
2. Open the shoe. The hook and loop fastener allows the Wound Care Shoe to be opened and closed quickly and easily.
3. The tried and tested 4-layer sole system is part of the delivery (Fig. 2) and can be processed individually (Fig. 6). All four insoles must always be in the shoe.
4. Position the foot with socks or bandage comfortably in the shoe.
5. Close the hook and loop straps so the foot is firmly secured.

Ulcer marking and cutting the soles:

1. Place a sterile film over the wound or pressure ulcer (Fig. 3) and draw a circle on the film around the affected area with a marker pen (Fig. 4).
2. The patient should now stand up on the insole in the shoe(s) so the ulcer marks the upper layer of the insole.

3. The next two layers (white, orange) are now punctured through the top layer (beige) using a pointed object, thus defining the cutout to be formed (Fig. 5).
4. Now draw the cutouts to be removed on the underside of the layers, as described in Fig. 7.
5. Cut out the marked area using a knife.

Sole cut:

- A No plantar open sores
- B Severity 0 and 1, Wagner scale, open wound, cutting angle around 30°
- C After wound healing (from B)
- D Severity of 2 and higher, Wagner scale, open wound
- E After wound healing (from D)
- F Transverse section of the relief of a pressure ulcer under the third metatarsal head

Plantar layer system:

	Material	Density	Thickness	Durometer*
	Plastazote®	soft	6.3 mm	26
	EVA	medium	6.3 mm	42
	EVA	firm	6.3 mm	50
	Insole foam	high	3.1 mm	29

* The insole durometers may vary within 2° to 3° of those listed.

Note: Modifications may only be carried out by specialised medical personnel. The period of use depends on the grade of the ulcer and should be coordinated with the treating doctor. For a physiologically correct gait, use the shoe in pairs if possible or with the TwinShoe® balance shoe for height balancing.

Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling / cleaning. The medical device is designed to be used for the usual duration of use prescribed by a doctor without showing a reduction in safety or performance. After the period of use prescribed by your doctor, check the function of the product together with your doctor or specialist dealer.

Care instructions

If the shoe is slightly soiled, you can wipe it with a damp cloth or clean it with a soft brush. Please do not use any aggressive cleaning agents (alcohol etc.). Do not put the shoes in the dryer or on a heat source, but leave them to dry at room temperature. The soles provided can be cleaned with lukewarm water and soap.

Disposal

 Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. Dispose of this product after use in accordance with applicable local laws and regulations.

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V-M Orthotics Ltd.
ALBO-Healthcare GmbH