WCS[®]





**WCS[®] / WCS[®]
Light Insoles***



- EN** Wound Care Shoe / * Four-Layer-Insole-System
DE Verband- und Wundtherapieschuh / * 4-Sohlen-System
ES Zapato para tratamiento de úlceras / * Sistema de cuatro plantillas
FR Chaussure pour le soin des plaies diabétiques / * Système de semelle à 4 couches
IT Calzatura ortopedica terapeutica / * Sistema a 4 solette
NL Wondbehandelingsschoen / * 4-zolen-systeem
DA Bandage- og sårsko / * 4-sålssystem
SV Sårbehandlingssko / * System med 4 sulor
PL But pooperacyjny i do terapii ran / * System 4-podeszwowy
SL Čevljev za nego ran / * Sistem s 4 podplati
HR Cipelna za liječenje rana / * Sustav s 4 uloška
PT Sapato para terapia de feridas / * Sistema de 4 palmilhas

DARCO

WCS[®], WCS[®] / WCS[®] Light Insoles

WCS® | Wound Care Shoe

WCS® / WCS® Light Insoles | Four-Layer-Insole-System

Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your doctor, specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation, etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This device assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product in direct contact with injured, irritated, damaged or broken skin. This device is not recommended for direct skin contact. You should use bandages, dressings or socks.

You must not drive a vehicle while you are wearing this product. Use extreme caution when climbing stairs (e.g. also stepping up onto curbs) and walking on uneven surfaces. Be particularly cautious when walking in wet, slippery and snowy conditions. Walk more slowly than usual.

Note: This device is allowed for single-patient use only. It does not need to be cleaned or reprocessed between uses.

Purpose

The Wound Care Shoe is mainly used for treating open and closed pressure sores and ulcerations (malum perforans). This enables specific offloading and pressure redistribution of the foot.

Indications / Contraindications

Indications

- For treatment of diabetic foot syndrome
- For healing of open wounds and pressure ulcers and ulcerations
- For pressure relief
- During the postoperative phase after toe and / or forefoot surgery
- Ulcerations up to Wagner grade 3
- Bony prominences

Contraindications

- Open fractures
- Ulcerations higher than Wagner grade 3

Design Features

The special sole construction of the Wound Care Shoe achieves optimum pressure redistribution. The multi-layer insole can be adapted individually to the patient to achieve selective offloading. The upper material is as seamless as possible to prevent pressing. It can be worn as an everyday shoe both indoors and outdoors.

Use

Modifications may only be carried out by specialised medical personnel. Determining the correct shoe size is very important for the functionality of the shoe.

1. Determine the correct shoe size from the insole length.
2. Open the shoe. The hook and loop fastener allows the Wound Care Shoe to be opened and closed quickly and easily.
3. The tried and tested 4-layer sole system is part of the delivery (Fig. 2) and can be processed individually (Fig. 6). All four insoles must always be in the shoe.
4. Position the foot with socks or bandage comfortably in the shoe.
5. Close the hook and loop straps so the foot is firmly secured.

Example for ulcer marking and cutting the soles:

1. Place a sterile film over the wound or pressure ulcer (Fig. 3) and draw a circle on the film around the affected area with a marker pen (Fig. 4).
2. The patient should now stand up on the insole in the shoe(s) so the ulcer marks the upper layer of the insole.
3. The next two layers (white, orange) are now punctured through the top layer (beige) using a pointed object, thus defining the cutout to be formed (Fig. 5: Example of a modification).

- Now draw the cutouts to be removed on the underside of the layers, as described in Fig. 7.
- Cut out the marked area using a knife.



Upper material modification (non-load-bearing sores):

- To remove the upper material (synthetic leather) without damaging the liner, pass a blunt object (i.e. paperclip) from one perforated hole to the next. Be careful not to puncture the liner (beige).
- Cut the synthetic leather using the blunt object as a guide.
- Insert a flat instrument between the synthetic leather and liner (beige) to continue the cut. Only cut within the given perforation holes.

Sole cut:

- A No plantar open sores
- B Severity 0 and 1, Wagner scale, open wound, cutting angle around 30°
- C After wound healing (from B)
- D Severity of 2 and higher, Wagner scale, open wound
- E After wound healing (from D)
- F Transverse section of the relief of a pressure ulcer under the third metatarsal head

Plantar layer system:

	Material	Density	Thickness	Durometer*
	Plastazote®	soft	6.3 mm	26
	EVA	medium	6.3 mm	42
	EVA	firm	6.3 mm	50
	Insole foam	high	3.1 mm	29

* The insole durometers may vary within 2° to 3° of those listed.

Note: Modifications may only be carried out by specialised medical personnel.

The period of use depends on the grade of the ulcer and should be coordinated with the treating doctor. For a physiologically correct gait, use the shoe in pairs if possible or with the TwinShoe® balance shoe for height balancing.

Combination with other products

- WCS® / WCS® Light Insoles
- TwinShoe®

Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling / cleaning. The medical device is designed to be used for the usual duration of use prescribed by a doctor without showing a reduction in safety or performance. After the period of use prescribed by your doctor, check the function of the product together with your doctor or specialist dealer.

Technical data

Sizes

XS (35 – 36), S (37 – 38), S / M (39 – 40), M (41 – 42), L (43 – 44),

XL (45 – 46), XXL (47 – 48)

Right and left version available.

Colour

Black

Material

Outsole: Rubber, EVA, bamboo

Insoles: Plastazote®, EVA, PU

Upper material: Polyester, Plastazote®, synthetic leather

Maintenance

The product does not require maintenance.

Care instructions

If the shoe is slightly soiled, you can wipe it with a damp cloth or clean it with a soft brush. Please do not use any aggressive cleaning agents (alcohol etc.). Do not put the shoes in the dryer or on a heat source, but leave them to dry at room temperature.

The soles provided can be cleaned with lukewarm water and soap.

Disposal

Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. Dispose of this product after use in accordance with applicable local laws and regulations.

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ALBO-Healthcare GmbH
DARCO UK Ltd