CHRISOFIX®  
— The flexi-fix healer —

## First aid wind-up splints

### Ref\_04

Aufrollbare Notfallschienen

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Férulas enrollables de primeros auxilios

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#### PRODUCT CHARACTERISTICS

Chrisofix® orthoses can ensure limb fixation according to diverse clinical requirements. Chrisofix® products are repeatedly adjustable, yet shape-retaining, are light and washable. They limit the movement of the limb only to the necessary degree, but when used with the corresponding strap they also immobilise the limb like a POP or thermoplastic splint. It is clinically proven that the application of Chrisofix® orthoses results in pain reduction and improvement of functioning capacity.

#### INDICATION

Recommended to use in a first aid situation for temporary fixation of an injured elbow, forearm, wrist, hand, knee, leg (cruris), ankle or foot. Recommended to compliment with bandaging.

#### APPLICATION

After unwinding the splint it needs to be adjusted to the shape required for the type of injury on the uninjured contralateral bodypart. It may be formed into a flat U-shaped bed/rill or a „sugar tong” shape. The splint needs to be complimented with a bandage. It may also be used as a supporting bandage e.g. in case of an ankle joint trauma.

#### MATERIAL AND PRODUCT CARE

Chrisofix® wind-up splints are made of a non-allergenic and phthalate-free sandwich composite consisting of polyethylene foam and aluminium. The layer that has direct contact with the skin is skin-friendly and soft, so there is no need for extra cushioning. The structure is repeatedly adjustable without heat or water, shape-retaining and X-ray transparent.

The orthosis is not a sterile device. It is washable at 30°C and should be left to air dry. Store at room temperature. None of our products contain any hazardous substances. Once used, please discard as household waste.

#### SAFETY PRECAUTIONS

The orthosis should be applied to clean, undamaged skin or following wound treatment.

Avoid an overly tight application or one that leaves mark on the skin.

The first application and major changes to the form of the orthosis should be performed by professionals.

Excessively frequent shaping / bending of the orthoses should be avoided as this may reduce the stability of the aluminium core. Chrisofix® orthoses are not suitable for heavy physical work. Patients should avoid any activity that may lead to the orthosis getting damaged, broken, wet, contaminated or having a major change in form. Failures resulting from misuse are not covered by the product warranty.

Do not use a broken or damaged orthosis.

Do not cut the product!

Avoid the combination with other fixation products on the same limb, except in case of medical instruction.

Chrisofix® products are produced according to the quality control system ISO:13485. However, should arise any complain regarding the product, please contact your local representative dealer. Please report any serious incident that has occurred in relation to the device to the competent authority of the Member State in which the user and/or patient is established, the manufacturer and the distributor.

Please keep the packaging and Instructions for Use until the end of the treatment.

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