
Thumb orthosis - Skier's thumb

Ref_23 - Right
Ref_24 - Left

- Daumenorthese - Skidaumen
- Orthèse du pouce - Pouce du skieur
- Tutore del pollice - Pollice dello sciatore
- Ortesis de pulgar - Pulgar del esquiador
- Hüvelykujj ortézis - UCL

PRODUCT CHARACTERISTICS

Chrisofix® orthoses can ensure limb fixation according to diverse clinical requirements. Chrisofix® products are repeatedly adjustable, yet shape-retaining, are light and washable. They limit the movement of the limb only to the necessary degree, but when used with the corresponding strap they also immobilise the limb like a POP or thermoplastic splint. It is clinically proven that the application of Chrisofix® orthoses results in pain reduction and improvement of functioning capacity.

INDICATION

Immobilisation of the I. metacarpophalangeal (I.MCP) joint in case of injuries, distortion, i.e. skier's thumb. Partial immobilization of the thumb saddle joint, recommended for day use. In addition, it is proposed for various thumb MCP injuries, distortions, rehabilitation following fractures, or as a complement to appropriate surgical procedures around the thumb MCP, at the discretion of the surgeon.

APPLICATION

Chrisofix® orthoses are supplied in a pre-shaped form ready for their application. Place the limb into the orthosis according to the illustration, and adjust the orthosis with gentle pressure. Finally, fix the limb inside the orthosis with the straps supplied.

MATERIAL AND PRODUCT CARE

Chrisofix® orthoses are made of a non-allergenic and phthalate-free sandwich composite consisting of polyethylene foam, polyamide, cotton, and corrugated and perforated aluminium. The layer that has direct contact with the skin is skin-friendly and soft, so there is no need for extra cushioning. The structure is breathable, repeatedly adjustable without heat or water, shape-retaining and X-ray transparent. The orthosis is not a sterile device. It is washable at 30 °C and should be left to air dry. Store at room temperature. None of our products contain any hazardous substances. Once used, please discard as household waste.

SAFETY PRECAUTIONS

The orthosis should be applied to clean, undamaged skin or following wound treatment. Avoid an overly tight application or one that leaves mark on the skin. The first application and major changes to the form of the orthosis should be performed by professionals. Excessively frequent shaping / bending of the orthoses should be avoided as this may reduce the stability of the aluminium core. Chrisofix® orthoses are not suitable for heavy physical work. Patients should avoid any activity that may lead to the orthosis getting damaged, broken, wet, contaminated or having a major change in form. Failures resulting from misuse are not covered by the product warranty. Do not use a broken or damaged orthosis. Do not cut the product! Avoid the combination with other fixation products on the same limb, except in case of medical instruction.

Chrisofix® products are produced according to the quality control system ISO:13485. However, should arise any complain regarding the product, please contact your local representative dealer. Please report any serious incident that has occurred in relation to the device to the competent authority of the Member State in which the user and/or patient is established, the manufacturer and the distributor.

Please keep the packaging and Instructions for Use until the end of the treatment.

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