Rib splint

Ref_60


Rippenschiene

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Férula costal

Bordasín

MAIN CHARACTERISTICS OF THE PRODUCT

Chrisofix® rib splint developed and patented worldwide for rib fractures and chest contusions is light and comfortable, and it does not restrict, but only supports adequately deep breathing. As a result, it is clinically proven to prevent subsequent complications and reduce the pain associated with each breath. Wearing the Chrisofix® rib splint helps restore respiratory function faster and may minimize the duration of a potential hospitalization.

Each Chrisofix® rib splint can only be applied once. It can be worn for approximately one week, and it is not necessary to change the splint during this period, it provides continuous analgesic effect while the patient is wearing it. One week is usually enough to overcome the painful period, but if longer treatment is required, another product should be applied to maintain its beneficial effects.

INDICATION

This product is indicated for chest contusions, rib fracture(s), serial or window fractures (to prevent paradoxical breathing), and following thoracic surgery.

APPLICATION

The fracture site(s) should be localized with X-ray prior to the application of the splint. The splint should be applied with the patient standing or potentially sitting, preferably not lying. Stand or sit straight as far as possible. Straight posture helps in the proper placement of the splint. Before applying the splint, hair and grease should be removed from the affected skin surface. Bending the product to the form of the chest should be performed on the other, non-injured side before applying on the injured region.

Align the splint to the chest so that the side marked with an arrow faces towards the chest and the arrow should be parallel to the ribs.

Peel off the protective paper layer from the surface of the rigid body of the splint.

Apply the splint on the area of the fractured rib(s) with the patient exhaling, so that the splint bridges the injured area.

After adhering the splint to the chest (but not before), peel off the protective paper layer from the bottom of the waterproof foil on all sides and, while inhaling, smooth the foil onto the skin surface of the chest.

Use the peel-off tab on the outer surface of the splint to remove the outer transparent protective foil. Thus, the adhered waterproof foil becomes permeable for air. Under no circumstances should this outer protective foil be left on the product and only remove it as the last step, otherwise the product becomes unusable.

CHARACTERISTICS OF THE MATERIAL AND HANDLING OF THE PRODUCT

The Chrisofix® rib splint is made of a phthalate- and allergen-free sandwich composite consisting of polyethylene foam, corrugated and perforated aluminium, skin-friendly adhesives, waterproof and air permeable foil. The materials used on the surface of the product are skin-friendly and non-allergenic. The structure of the product, together with the waterproof foil, is breathable, form-retaining and radiolucent, so it does not cover the injury when taking follow-up X-rays, so the healing of the injury can be monitored on the X-ray images even with the patient wearing the splint. Its waterproof cover allows the injured person to take a shower (but not to go in a bathtub or a pool), making everyday personal hygiene more comfortable. The product is not sterile. It does not contain hazardous substances; it can be disposed of with household waste after use. Recommended storage: 10 -27 °C and 40-60 % humidity.

WARNING

The orthosis can be applied on clean, undamaged skin or after wound care. It should be applied with the assistance of a specialist.

Activities that cause damage, breakage, significant deformation or contamination of the product must be avoided.

Do not use a product that is broken or whose structure is damaged.

The product must not be cut!

In rare cases, itching may occur; in this case, it is advisable to remove the foil covering the splint and to enhance the fixation of the splint to the skin with a non-allergenic adhesive plaster. Once the cover foil has been removed, showering should only take place after a temporarily applied waterproof cover. Avoid using multiple fixation products at the same time, unless recommended by your physician.

Chrisofix® products are manufactured according to the strict rules of the ISO 13485 quality management system. Should nonetheless any complaint arise regarding the product, you can report it to your local distributor. Please report any serious incidents that has occurred in relation to the device to the competent authority of the Member State in which the user and/or patient is established, the manufacturer and the distributor.

Please keep the packaging and Instructions for Use until the end of the treatment.

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