# DARCO Posterior Splint Fix

EN Foot-positioning orthosis  
DE Fußlagerungsorthese  
ES Ortesis estabilizadora del pie  
FR Orthèse de positionnement du pied  
IT Ortesi da riposo  
NL Voetimmobilisatiebrace  
DA Fodlejringsortose  
SV Viloortos  
PL Orteza unieruchamiająca stopę  
SL Položajna ortoza za stopalo  
HR Imobilizacijska ortoza za stopalo  
PT Órtese de posicionamento do pé

**DARCO**

1.

2.

3.

4.

### DARCO Posterior Splint Fix | Foot-positioning orthosis

#### Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This aid assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product with injured, irritated, damaged or intact skin. You should use bandages or socks, for example.

This aid is intended for single-patient use only.

## Purpose

Foot-positioning orthosis for immobilization in a defined position. For use as a post- and preoperative positioning orthosis, for immobilisation following accidents, post-traumatic immobilisation and instead of plaster.

## Indications / Contraindications

### Indications

- For general use in conservative and postoperative therapy of lower leg, ankle and foot injuries:
- > Postoperative stabilisation
  - > Post-traumatic stabilisation / immobilisation
  - > Immobilisation following an accident
  - > Short-term stabilisation (ligament rupture, luxation, fractures, wounds, strains, sprains, etc.)
  - > Preoperative stabilisation

110° models

- > Achilles tendon rupture

### Contraindications

- > Calf muscle shortening
- > Marked Haglund / posterolateral exostosis
- > Impaired sensation (e.g. polyneuropathy), vascular disorders
- > Material incompatibility / skin changes (e.g. swelling, redness) due to allergic reactions, inflammation or injury.

## Design features

Stable, shelled plastic orthosis intended for positioning in a defined position with comfortable liner and flexibly adjustable fixing straps.

## Combination with other products

No specifications

## Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling / cleaning.

**Putting on the positioning orthosis:**

- > First treat the surgical wound or foot injury in the usual manner with appropriate bandaging material.
- > Select the appropriate length of the orthosis. The patient's foot and lower leg length is the determining factor for this.

**Note**

On the lower leg, it is important that the orthosis does not extend too far towards the back of the knee in order to avoid pressure in this area when flexing the knee joint. It should also be ensured that the head of the fibula is exposed in order to avoid pressure damage to the peroneal nerve. Do not wear directly on open wounds or irritated skin.

- > Position the foot in the orthosis. To do this, the heel must be correctly positioned in the heel recess. (photo 1)
- > Place the wide strap around the calf area of the orthosis and close it in a circle. (photo 2)

Ensure that the hook strap does not rest on the skin.

- > Wind a narrow strap around the patient's forefoot. For more stability, fix the second narrow strap in the desired position. (photo 3) The strap material adheres to the hook and loop straps of the shell. (photo 4) Ensure that no hook strap rests on the skin.

**Alternative use:**

- > Secure the orthosis using a conventional gauze bandage or cohesive bandage. Integrated hook and loop straps on the outer shell allow ideal adhesion of the gauze bandage. Avoid winding too tightly.
- > The bandage material should ideally be approx. 10 cm wide. Begin winding the bandage at the back of the foot and guide the bandage in a distal (away from body) to a proximal direction (towards body). In doing so, check that the orthosis is correctly positioned.

## Technical data

<b>Sizes:</b>	Foot length of the orthosis	Calf section of the orthosis
S 90° Plus/S 110°	265 mm	340 mm
M 90° Plus/M 110°	295 mm	370 mm
L 90° Plus/L 110°	325 mm	400 mm

Can be used universally on the right or left. The 110° models have been developed for Achilles tendon pathologies in a 20° club foot position.

**Colour:** White-blue

**Material**

<b>Outer shell:</b>	Polycarbonate (thermoplastic deformable)
<b>Inner lining:</b>	PE foam
<b>Hook and loop straps:</b>	Nylon
<b>Fixation straps:</b>	Viscose, elastane, polyamide

## Maintenance

The product does not require maintenance.

## Cleaning instructions

- > Wipe disinfection of the outer shell and inner liner with a surface disinfectant
- > Wash the straps at 30°C. Close the straps for washing.
- > Do not use any substances that contain oil

## Disposal

Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. The product can be disposed of with household waste.

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