New Yogipis Plus English ver.	
Instructions for use	New Yogipis Plus is intended for use in most common orthopedic splint treatments for outpatient, rehabilitation, emergency, before and after surgery, replacement of casting tape, and injuries such as ligament injuries, muscle rupture, sprain and strain damage.
Contraindications	No contraindication known at present time.
Application Instructions	The splint can be formed with the help of a heat source. To do this, proceed as following: 1. Apply heat by using any heating device, such as Woosam Cardboard Oven, heat gun, hair dryer & etc. - If you only want to mold a specific area of the product, it is recommended to heat only this area and not the whole product. 2. After 2-3 minutes of heating, the product becomes moldable and cuttable. 3. If needed, you can now cut the product to size with an appropriate scissors. In any case, note that the edges may be sharp after cutting and sufficient padding must therefore be attached to the edges. 4. Place the product on the relevant part of the patient's body. Make sure that the product has cooled down sufficiently after heating and before putting it on the patient. Otherwise there is a risk of burns. 5. Shape the product to the patient's anatomy. 6. Fix the still flexible splint for example with bandages so that the splint keeps its shape. 7. After 2-5 minutes the product stops being moldable and becomes hard. 8. Remove any bandages attached to fix the splint. 9. Fix the splint by closing the straps with the Velcro fasteners.
Important Notes	1. To avoid overheating, the workplace must be well-ventilated. 2. Keep away from open flame. 3. If necessary, especially in the presence of sharp edges, use additional padding. 4. This product is intended for multiple use on a single patient only. 5. Only by choosing the correct size and variant, as well as the correct application, an optimal function can be achieved. 6. The product should no longer be used if it is damaged or heavily worn, or if existing joints and fasteners no longer function properly. 7. Make sure that existing Velcro fasteners do not cause damage to clothes or other fabric parts. 8. Do not put on the product too tightly, as this may affect the blood circulation. 9. If symptoms persist or worsen, please consult your medical doctor. 10. If you are unclear or have questions about the product and its handling, please consult your doctor or qualified personnel. 11. Do not make any changes to the product. 12. Simultaneous use with other products only after consultation with your doctor. 13. Any serious incident (death, serious deterioration in health, or a threat to public health) occurring in connection with this medical device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is located.
Material composition	Hard outer shell: Thermoplastic Padding: Air mesh fabric Straps: Polyamide, Polyester, Nylon
Storage Instruction	Store the product inside in a dry place and protect it from heat, direct sunlight and moisture. Storage at room temperature (15-30 degrees Celsius) is recommended.
Disposal note	Disposal via household waste is possible. It is recommended to follow local regulatory laws for proper disposal
Shelf Life	3years (No specific long-term stability testing has been conducted. Properly stored products may remain usable beyond the stated shelf life.
  	Orthopartner AG Lenzburgerstrasse 2 Gebäude 10 5702 Niederlenz Switzerland
	AF PHARMA SERVICE EUROPE SL (SRN ES-AR-000018076) Muntaner 281, 08021, Barcelona, España info@alpharmaservice.com, T. +34690846446
	Woosam Medical Co., Ltd. (SRN: KR-MF-000045208) 89-25, Tongil-ro642beon-gil, Wollong-myeon, Paju-si, Gyeonggi-do, 10945, Rep. of Korea
	Woosam Medical Co., Ltd. (SRN: KR-MF-000045208) 89-25, Tongil-ro642beon-gil, Wollong-myeon, Paju-si, Gyeonggi-do, 10945, Rep. of Korea
     	
Version: V02/2026/06	

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