Hygia Longuette English ver.	
Intended Use	Longuettes are designed for use either independently as individual splints or together with Woosam's rigid or semi-rigid casting tapes for added support.
Contraindications	Currently unknown.
Application Instructions	1. Begin by covering the affected area with one or two layers of Stockinette, then wrap it with Padding for skin protection. 2. Always wear appropriate gloves when handling the Longuette to prevent resin transfer to your skin. Uncured resin can stick to the skin or clothing, so avoid any direct contact with the patient's skin until the material has fully cured or is properly covered. 3. If applying the Longuette alongside casting tape, start by applying the initial layer of casting tape first. 4. Choose the suitable size of the Longuette. Before opening, gently squeeze the sealed pouch to ensure the material inside is soft and pliable. Do not use any splint that feels partially or completely hardened. After opening, the splint can be trimmed to the desired shape and size with bandage scissors if customization is needed. 5. Submerge the Longuette in water at room temperature (21-24°C). Once soaked, remove it and gently squeeze out the excess water. Repeated squeezing while immersed can speed up the setting time without the need to increase the water temperature. During curing, the Longuette will generate heat due to an exothermic reaction. When used with circular casting tape, the Longuette may be applied dry and secured with the subsequent layer of casting tape, either wet or dry. Additionally, applying a damp elastic bandage over the synthetic casting layers will also initiate polymerization. 6. Place the Longuette onto the affected limb and secure it by wrapping additional layers of casting tape. 7. Shape and mold the splint immediately after application. Hold the splint firmly in areas requiring detailed molding for approximately 30 seconds to accelerate hardening. - Keep the limb in the correct position until the material has fully set, usually within three to four minutes. 8. Note that a splint alone is not meant for weight-bearing use. However, when the Longuette is combined with Soft Cast material for reinforcement, it will provide sufficient weight-bearing support after 20 to 30 minutes.
Care & Precaution	The uncured resin adheres firmly to skin and to clothing. Swabbing with alcohol or warm soapy water will be helpful to remove uncured resin immediately when touched to skin or clothing. Shower is not recommended with the splint on. Follow a doctor's instructions for taking shower.
Reporting Obligation	Any serious incident (death, serious deterioration in health, or a threat to public health) occurring in connection with this medical device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.
Material composition	Non woven fabric Knitted polyester fabric Polyurethane resin
Storage Instruction	Store in a dry and cool condition at temperature between 15-30°C. Extremely high or low temperature may reduce shelf life. Care should be taken to avoid puncturing the pouch as this will cause the tape to harden prematurely in the package. Refer to the expiry date shown on products and labels attached.
Disposal note	Please respect the local laws and regulation for waste disposal.
Shelf Life	3years
	Orthopartner AG Lenzburgerstrasse 2 Gebäude 10 5702 Niederlenz Switzerland
	AF PHARMA SERVICE EUROPE SL (SRN ES-AR-000018076) Muntaner 281, 08021, Barcelona, España info@afpharmaservice.com, T. +34690846446
	Woosam Medical Co., Ltd. (SRN: KR-MF-000045208) 89-25, Tongil-ro642beon-gil, Wollong-myeon, Paju-si, Gyeonggi-do, 10945, Rep. of Korea
     Version: V01/2025/07	

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