


Oval-8[®] Finger Splints

- Stabilizes or immobilizes IP finger joints
- Simple to fit and wear for acute or long term use
- Open 3-point pressure design allows for treatment of multiple conditions
- Allows patients to maintain full hand function—can be worn in water
- Can be adjusted in seconds for a custom fit using a heat gun



Oval-8[®] Finger Splints - Indications



Mallet Finger



Trigger Finger



Trigger Thumb



Swan Neck Deformity



Boutonniere Deformity



Lateral Deviation



Immobilization



Oval-8[®] Finger Splint New Look



Oval-8[®] Finger Splints

The Best are Even Better

Oval-8 splints are now easier on your eyes and easier on your patient's finger



The Oval-8 name and size are now **easier to read**. Hurray!



The new **+** sign marks the end that fits looser when put on the finger first



The seamless bands are **smoother & more comfortable** on injured fingers

Oval-8[®] Finger Splint – Buying Options



Oval-8 Kits

- *Oval-8 Kit*
- *Oval-8 Pediatric Kit*



Oval-8 Sizing Set

Great for small clinics that do not dispense splints



Single & Mixed Size Packages

- *Packages of 1*
- *Packages of 5*
- *Combo Packages*
- *Graduated Sets*